




FCRL4 Antibody

Product Code	CSB-PA836277LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PJ5
Immunogen	Recombinant Human Fc receptor-like protein 4 protein (20-387AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fc receptor-like protein 4 (FcR-like protein 4) (FcRL4) (Fc receptor homolog 4) (FcRH4) (IFGP family protein 2) (hIFGP2) (Immune receptor translocation-associated protein 1) (CD antigen CD307d), FCRL4, FCRH4 IFGP2 IRTA1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FCRL4

Image



Western blot
 All lanes: FCRL4 antibody at 6µg/ml + 293T whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 58, 57 kDa
 Observed band size: 58 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.