



FCRL4 Antibody, Biotin conjugated

Product Code	CSB-PA836277LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PJ5
Immunogen	Recombinant Human Fc receptor-like protein 4 protein (20-387AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fc receptor-like protein 4 (FcR-like protein 4) (FcRL4) (Fc receptor homolog 4) (FcRH4) (IFGP family protein 2) (hIFGP2) (Immune receptor translocation-associated protein 1) (CD antigen CD307d), FCRL4, FCRH4 IFGP2 IRTA1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FCRL4