



FPR1 Antibody, Biotin conjugated

Product Code	CSB-PA07129D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21462
Immunogen	Recombinant Human fMet-Leu-Phe receptor protein (306-350AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	fMet-Leu-Phe receptor (fMLP receptor) (N-formyl peptide receptor) (FPR) (N-formylpeptide chemoattractant receptor), FPR1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FPR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.