



FSCB Antibody, HRP conjugated

Product Code	CSB-PA708110LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5H9T9
Immunogen	Recombinant Human Fibrous sheath CABYR-binding protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibrous sheath CABYR-binding protein, FSCB, C14orf155
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FSCB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.