



GEMIN8 Antibody, HRP conjugated

Product Code	CSB-PA873675LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWZ8
Immunogen	Recombinant Human Gem-associated protein 8 protein (1-242AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gem-associated protein 8 (Gemin-8) (Protein FAM51A1), GEMIN8, FAM51A1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GEMIN8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.