



GMEB1 Antibody, Biotin conjugated

Product Code	CSB-PA896751LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y692
Immunogen	Recombinant Human Glucocorticoid modulatory element-binding protein 1 protein (2-245AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucocorticoid modulatory element-binding protein 1 (GMEB-1) (DNA-binding protein p96PIF) (Parvovirus initiation factor p96) (PIF p96), GMEB1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GMEB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.