





GMPPB Antibody

Product Code	CSB-PA896740LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5P6
Immunogen	Recombinant Human Mannose-1-phosphate guanyltransferase beta protein (61-360AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mannose-1-phosphate guanyltransferase beta (EC 2.7.7.13) (GDP-mannose pyrophosphorylase B) (GTP-mannose-1-phosphate guanylyltransferase beta), GMPPB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GMPPB

Image



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA896740LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA896740LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.