



GPRC5B Antibody, HRP conjugated

Product Code	CSB-PA868367LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZH0
Immunogen	Recombinant Human G-protein coupled receptor family C group 5 member B protein (293-403AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	G-protein coupled receptor family C group 5 member B (A-69G12.1) (Retinoic acid-induced gene 2 protein) (RAIG-2), GPRC5B, RAIG2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GPRC5B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.