





GSDMB Antibody

Product Code	CSB-PA823445LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TAX9
Immunogen	Recombinant Human Gasdermin-B protein (63-310AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gasdermin-B (Gasdermin-like protein), GSDMB, GSDML
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GSDMB

Image



Western blot
All lanes: GSDMB antibody at 8µg/ml + Colo320 whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 47, 46, 48, 19 kDa
Observed band size: 47 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA823445LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.