



HDHD3 Antibody, HRP conjugated

Product Code	CSB-PA883592LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BSH5
Immunogen	Recombinant Human Haloacid dehalogenase-like hydrolase domain-containing protein 3 protein (1-251AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Haloacid dehalogenase-like hydrolase domain-containing protein 3, HDHD3, C9orf158
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HDHD3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.