




HMGN3 Antibody

Product Code	CSB-PA623009LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15651
Immunogen	Recombinant Human High mobility group nucleosome-binding domain-containing protein 3 protein (1-77AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	High mobility group nucleosome-binding domain-containing protein 3 (Thyroid receptor-interacting protein 7) (TR-interacting protein 7) (TRIP-7), HMGN3, TRIP7
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HMGN3

Image


Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA623009LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.