



ICA1L Antibody, Biotin conjugated

Product Code	CSB-PA847691LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NDH6
Immunogen	Recombinant Human Islet cell autoantigen 1-like protein (1-187AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Islet cell autoantigen 1-like protein (Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 14 protein) (Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 15 protein), ICA1L, ALS2CR14 ALS2CR15
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ICA1L