



IER2 Antibody

Product Code	CSB-PA01729A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BTL4
Immunogen	Recombinant Human Immediate early response gene 2 protein (1-223AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Immediate early response gene 2 protein (Protein ETR101), IER2, ETR101 PIP92
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	IER2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.