





KATNAL2 Antibody

Product Code	CSB-PA809022LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IYT4
Immunogen	Recombinant Human Katanin p60 ATPase-containing subunit A-like 2 protein (18-307AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Katanin p60 ATPase-containing subunit A-like 2 (Katanin p60 subunit A-like 2) (EC 3.6.4.3) (p60 katanin-like 2), KATNAL2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KATNAL2

Image



Western blot
 All lanes: KATNAL2 antibody at 1µg/ml + U251 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 62, 53 kDa
 Observed band size: 62 kDa



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA809022LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.