



KCNRG Antibody, Biotin conjugated

Product Code	CSB-PA847638LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N5I3
Immunogen	Recombinant Human Potassium channel regulatory protein (1-229AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Potassium channel regulatory protein (Potassium channel regulator) (Protein CLLD4), KCNRG, CLLD4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNRG