



KDM8 Antibody, Biotin conjugated

Product Code	CSB-PA839789LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N371
Immunogen	Recombinant Human JmjC domain-containing protein 5 protein (151-380AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	JmjC domain-containing protein 5 (EC 1.14.11.27) (EC 3.4.-.-) (Jumonji domain-containing protein 5) (Lysine-specific demethylase 8), KDM8, JMJD5
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	KDM8