



KDM8 Antibody, Biotin conjugated

Product Code	CSB-PA839789LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N371
Immunogen	Recombinant Human JmjC domain-containing protein 5 protein (151-380AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	JmjC domain-containing protein 5 (EC 1.14.11.27) (EC 3.4.-.-) (Jumonji domain-containing protein 5) (Lysine-specific demethylase 8), KDM8, JMJD5
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	KDM8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.