




KHDC4 Antibody

Product Code	CSB-PA770364LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z7F0
Immunogen	Recombinant Human KH homology domain-containing protein 4 protein (1-241AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	KH homology domain-containing protein 4 (Brings lots of money 7) (Pre-mRNA splicing factor protein KHDC4), KHDC4, BLOM7 KIAA0907 SNORA80EHG
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	KHDC4

Image



Western blot


All lanes: KHDC4 antibody at 5µg/ml + Mouse liver tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 65, 60, 41, 25 kDa

Observed band size: 65 kDa



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA770364LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.