



KLHDC8A Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA816911LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8IYD2 |
| Immunogen | Recombinant Human Kelch domain-containing protein 8A protein (1-350AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Kelch domain-containing protein 8A (Substitute for delta-EGFR expression 1) (S-delta-E1), KLHDC8A |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | KLHDC8A |