



Klk1b22 Antibody

Product Code	CSB-PA325391LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15948
Immunogen	Recombinant Mouse Kallikrein 1-related peptidase b22 protein (25-259AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kallikrein 1-related peptidase b22 (EC 3.4.21.35) (Beta-NGF-endopeptidase) (Epidermal growth factor-binding protein type A) (EGF-BP A) (Glandular kallikrein K22) (mGK-22) (Nerve growth factor beta chain endopeptidase) (Tissue kallikrein 22), Klk1b22, Klk-22 Klk22
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Klk1b22
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.