



# Klk1b22 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA325391LB01MO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P15948
<b>Immunogen</b>	Recombinant Mouse Kallikrein 1-related peptidase b22 protein (25-259AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Kallikrein 1-related peptidase b22 (EC 3.4.21.35) (Beta-NGF-endopeptidase) (Epidermal growth factor-binding protein type A) (EGF-BP A) (Glandular kallikrein K22) (mGK-22) (Nerve growth factor beta chain endopeptidase) (Tissue kallikrein 22), Klk1b22, Klk-22 Klk22
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Research Area</b>	Others
<b>Target Names</b>	Klk1b22