



L3HYPDH Antibody, HRP conjugated

Product Code	CSB-PA846613LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EM0
Immunogen	Recombinant Human Trans-3-hydroxy-L-proline dehydratase protein (13-312AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trans-3-hydroxy-L-proline dehydratase (EC 4.2.1.77) (Trans-L-3-hydroxyproline dehydratase), L3HYPDH, C14orf149
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	L3HYPDH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.