




LGALS12 Antibody

Product Code	CSB-PA846606LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96DT0
Immunogen	Recombinant Human Galectin-12 protein (1-336AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Galectin-12 (Gal-12) (Galectin-related inhibitor of proliferation), LGALS12, GRIP1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LGALS12

Image



Western Blot

Positive WB detected in: A549 whole cell lysate, Rat kidney tissue


All lanes: LGALS12 antibody at 3µg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 38, 32, 37, 30, 36 kDa

Observed band size: 38 kDa



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA846606LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA846606LA01HU at dilution of 1:100