




MACROD1 Antibody

Product Code	CSB-PA860244LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BQ69
Immunogen	Recombinant Human O-acetyl-ADP-ribose deacetylase MACROD1 protein (83-325AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	O-acetyl-ADP-ribose deacetylase MACROD1 (EC 3.2.2.-) (EC 3.5.1.-) (MACRO domain-containing protein 1) (Protein LRP16) ([Protein ADP-ribosylglutamate] hydrolase), MACROD1, LRP16
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MACROD1

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA860244LA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.