





MAGEA9 Antibody

Product Code	CSB-PA013335LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43362
Immunogen	Recombinant Human Melanoma-associated antigen 9 protein (1-315AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Melanoma-associated antigen 9 (Cancer/testis antigen 1.9) (CT1.9) (MAGE-9 antigen), MAGEA9; MAGEA9B, MAGE9 MAGEA9A;
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	MAGEA9

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA013335LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA013335LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.