



MED18 Antibody, FITC conjugated

Product Code	CSB-PA871577LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BUE0
Immunogen	Recombinant Human Mediator of RNA polymerase II transcription subunit 18 protein (35-208AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mediator of RNA polymerase II transcription subunit 18 (Mediator complex subunit 18) (p28b), MED18
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MED18