



MIF4GD Antibody

Product Code	CSB-PA013827LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A9UHW6
Immunogen	Recombinant Human MIF4G domain-containing protein (1-133AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MIF4G domain-containing protein (SLBP-interacting protein 1) (hSLIP1), MIF4GD, SLIP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MIF4GD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.