



MRPL33 Antibody, Biotin conjugated

Product Code	CSB-PA014839LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75394
Immunogen	Recombinant Human 39S ribosomal protein L33, mitochondrial protein (1-65AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	39S ribosomal protein L33, mitochondrial (L33mt) (MRP-L33) (Mitochondrial large ribosomal subunit protein bL33m), MRPL33, C2orf1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MRPL33
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.