



MRPL33 Antibody, HRP conjugated

Product Code	CSB-PA014839LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75394
Immunogen	Recombinant Human 39S ribosomal protein L33, mitochondrial protein (1-65AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	39S ribosomal protein L33, mitochondrial (L33mt) (MRP-L33) (Mitochondrial large ribosomal subunit protein bL33m), MRPL33, C2orf1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MRPL33