



# MRPS18C Antibody

<b>Product Code</b>	CSB-PA014893LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9Y3D5
<b>Immunogen</b>	Recombinant Human 28S ribosomal protein S18c, mitochondrial protein (1-142AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	28S ribosomal protein S18c, mitochondrial (MRP-S18-c) (Mrps18-c) (S18mt-c) (28S ribosomal protein S18-1, mitochondrial) (MRP-S18-1) (Mitochondrial small ribosomal subunit protein bS18c) (Mitochondrial small ribosomal subunit protein bS18m), MRPS18C
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Tags & Cell Markers
<b>Target Names</b>	MRPS18C
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.