



MRPS18C Antibody, Biotin conjugated

Product Code	CSB-PA014893LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3D5
Immunogen	Recombinant Human 28S ribosomal protein S18c, mitochondrial protein (1-142AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	28S ribosomal protein S18c, mitochondrial (MRP-S18-c) (Mrps18-c) (S18mt-c) (28S ribosomal protein S18-1, mitochondrial) (MRP-S18-1) (Mitochondrial small ribosomal subunit protein bS18c) (Mitochondrial small ribosomal subunit protein bS18m), MRPS18C
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MRPS18C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.