



MS4A3 Antibody, Biotin conjugated

Product Code	CSB-PA856947LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96HJ5
Immunogen	Recombinant Human Membrane-spanning 4-domains subfamily A member 3 protein (1-49AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Membrane-spanning 4-domains subfamily A member 3 (CD20 antigen-like protein) (Hematopoietic-specific transmembrane protein 4) (HTm4), MS4A3, CD20L HTM4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MS4A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.