



NAPA Antibody, FITC conjugated

Product Code	CSB-PA015447LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54920
Immunogen	Recombinant Human Alpha-soluble NSF attachment protein (1-295AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-soluble NSF attachment protein (SNAP-alpha) (N-ethylmaleimide-sensitive factor attachment protein alpha), NAPA, SNAPA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NAPA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.