



NAT2 Antibody, Biotin conjugated

Product Code	CSB-PA015471LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11245
Immunogen	Recombinant Human Arylamine N-acetyltransferase 2 protein (1-290AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Arylamine N-acetyltransferase 2 (EC 2.3.1.5) (Arylamide acetylase 2) (N-acetyltransferase type 2) (NAT-2) (Polymorphic arylamine N-acetyltransferase) (PNAT), NAT2, AAC2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NAT2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.