



# Nkx2-2 Antibody, FITC conjugated

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA015842LC01MO  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                     |
| <b>Uniprot No.</b>         | P42586  |
| <b>Immunogen</b>           | Recombinant Mouse Homeobox protein Nkx-2.2 protein (1-273AA)                      |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Mouse   |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | FITC  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4  |
| <b>Purification Method</b> | >95%, Protein G purified  |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Alias</b>               | Homeobox protein Nkx-2.2 (Homeobox protein NK-2 homolog B), Nkx2-2, Nkx-2.2 Nkx2b |
| <b>Immunogen Species</b>   | Mus musculus (Mouse)  |
| <b>Research Area</b>       | Neuroscience  |
| <b>Target Names</b>        | Nkx2-2  |