



NME7 Antibody, Biotin conjugated

Product Code	CSB-PA897289LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5B8
Immunogen	Recombinant Human Nucleoside diphosphate kinase 7 protein (1-150AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoside diphosphate kinase 7 (NDK 7) (NDP kinase 7) (EC 2.7.4.6) (nm23-H7), NME7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NME7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.