



NOL6 Antibody, Biotin conjugated

Product Code	CSB-PA872506LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H6R4
Immunogen	Recombinant Human Nucleolar protein 6 protein (951-1146AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleolar protein 6 (Nucleolar RNA-associated protein) (Nrap), NOL6
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NOL6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.