



NRG2 Antibody, Biotin conjugated

Product Code	CSB-PA016078LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14511
Immunogen	Recombinant Human Pro-neuregulin-2, membrane-bound isoform protein (112-405AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pro-neuregulin-2, membrane-bound isoform (Pro-NRG2) [Cleaved into: Neuregulin-2 (NRG-2) (Divergent of neuregulin-1) (DON-1) (Neural- and thymus-derived activator for ERBB kinases) (NTAK)], NRG2, NTAK
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NRG2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.