



NSUN4 Antibody, HRP conjugated

Product Code	CSB-PA850270LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96CB9
Immunogen	Recombinant Human 5-methylcytosine rRNA methyltransferase NSUN4 protein (26-284AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	5-methylcytosine rRNA methyltransferase NSUN4 (EC 2.1.1.-) (NOL1/NOP2/Sun domain family member 4), NSUN4
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NSUN4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.