



NUDT9 Antibody, FITC conjugated

Product Code	CSB-PA880117LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BW91
Immunogen	Recombinant Human ADP-ribose pyrophosphatase, mitochondrial protein (47-350AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ADP-ribose pyrophosphatase, mitochondrial (EC 3.6.1.13) (ADP-ribose diphosphatase) (ADP-ribose phosphohydrolase) (Adenosine diphosphoribose pyrophosphatase) (ADPR-PPase) (Nucleoside diphosphate-linked moiety X motif 9) (Nudix motif 9), NUDT9, NUDT10
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NUDT9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.