



NUP37 Antibody, Biotin conjugated

Product Code	CSB-PA847704LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFH4
Immunogen	Recombinant Human Nucleoporin Nup37 protein (1-326AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoporin Nup37 (p37) (Nup107-160 subcomplex subunit Nup37), NUP37
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NUP37
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.