




OSBPL2 Antibody

Product Code	CSB-PA884442LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1P3
Immunogen	Recombinant Human Oxysterol-binding protein-related protein 2 protein (251-480AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oxysterol-binding protein-related protein 2 (ORP-2) (OSBP-related protein 2), OSBPL2, KIAA0772 ORP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	OSBPL2

Image



Western blot

All lanes: OSBPL2 antibody at 1µg/ml

Lane 1: HepG2 whole cell lysate

Lane 2: Mouse heart tissue


Lane 3: Hela whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 56, 54 kDa

Observed band size: 56 kDa



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA884442LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.