



PBLD Antibody, HRP conjugated

Product Code	CSB-PA017501LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30039
Immunogen	Recombinant Human Phenazine biosynthesis-like domain-containing protein (11-288AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phenazine biosynthesis-like domain-containing protein (EC 5.1.-.-) (MAWD-binding protein) (MAWDBP) (Unknown protein 32 from 2D-page of liver tissue), PBLD, MAWBP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PBLD