



PCBP4 Antibody, HRP conjugated

Product Code	CSB-PA017521LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57723
Immunogen	Recombinant Human Poly(rC)-binding protein 4 protein (1-150AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Poly(rC)-binding protein 4 (Alpha-CP4), PCBP4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PCBP4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.