



PHF10 Antibody, Biotin conjugated

Product Code	CSB-PA017898LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUB8
Immunogen	Recombinant Human PHD finger protein 10 protein (151-360AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PHD finger protein 10 (BRG1-associated factor 45a) (BAF45a) (XAP135), PHF10, BAF45A
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PHF10