





PITPNC1 Antibody

Product Code	CSB-PA892330LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKF7
Immunogen	Recombinant Human Cytoplasmic phosphatidylinositol transfer protein 1 protein (1-268AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytoplasmic phosphatidylinositol transfer protein 1 (Mammalian rdgB homolog beta) (M-rdgB beta) (MrdgBbeta) (Retinal degeneration B homolog beta) (RdgBbeta), PITPNC1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PITPNC1

Image



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA892330LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA892330LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.