



PKD1L2 Antibody, FITC conjugated

Product Code	CSB-PA773589LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z442
Immunogen	Recombinant Human Polycystic kidney disease protein 1-like 2 protein (1-306AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polycystic kidney disease protein 1-like 2 (PC1-like 2 protein) (Polycystin-1L2), PKD1L2, KIAA1879 PC1L2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PKD1L2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.