





PPM1M Antibody

Product Code	CSB-PA842717LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96MI6
Immunogen	Recombinant Human Protein phosphatase 1M protein (1-247AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein phosphatase 1M (EC 3.1.3.16) (Protein phosphatase 2C isoform eta) (PP2C-eta) (PP2CE), PPM1M, PPM1E
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPM1M

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA842717LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA842717LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.