



PPM1M Antibody, Biotin conjugated

Product Code	CSB-PA842717LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96MI6
Immunogen	Recombinant Human Protein phosphatase 1M protein (1-247AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein phosphatase 1M (EC 3.1.3.16) (Protein phosphatase 2C isoform eta) (PP2C-eta) (PP2CE), PPM1M, PPM1E
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPM1M
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.