



PTPA Antibody, FITC conjugated

Product Code	CSB-PA620876LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15257
Immunogen	Recombinant Human Serine/threonine-protein phosphatase 2A activator protein (1-250AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein phosphatase 2A activator (EC 5.2.1.8) (PP2A, subunit B', PR53 isoform) (Phosphotyrosyl phosphatase activator) (PTPA) (Serine/threonine-protein phosphatase 2A regulatory subunit 4) (Serine/threonine-protein phosphatase 2A regulatory subunit B'), PTPA, PPP2R4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PTPA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.